# Virtual® 380

VPS Impression Materials

Instructions for Use

Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Brugsanvisning

- Vinylpolysiloxane impression material
- Vinylpolysiloxan-Abformmasse (additionsvernetzend)
- Matériau d'empreinte à base de vinyle polysiloxane
- Materiale per impronte al polivinilsilossano (silicone d'addizione)
- Material de impresión de Vinilpolisiloxano
- Vinilpolisiloxano para moldagem
- Vinylpolysiloxan avtrycksmaterial
- Vinylpolysiloxan aftryksmateriale

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Bruksanvisning

Productinformatie

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For dental use only!

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist.

600156/0811/WE3

ivoclar  
vivadent®  
clinical

# Instructions for Use

## Description

Virtual impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

## Colors

See table "Technical Data"

## Working times and setting times

Material	Speed	Total Work Time [min:sec]	Mouth Set Time* [min:sec]
Monophase [tray/wash]	Fast Set	2:15	3:30
Heavy Body [tray]	Regular Set	2:05	4:30
	Fast Set	1:15	2:30

– Set times may deviate from the indicated values by (+/-) 15 seconds.

\* Minimum time the impression material should remain in the mouth before removal.

## Technical Data

	Monophase	Heavy Body
Color	Blue	Blue
Classification ISO 4823 / ADA Spec. No. 19	Type 2 Medium-bodied consistency	Type 1 Heavy-bodied consistency
Mixing Ratio [Base : Catalyst]	5:1	5:1
Linear Dimensional Change [24h]	< 0.20 %	< 0.20 %
Strain in Compression [min. – max.]	3 – 5 %	3 – 5 %
Recovery from Deformation	> 99.5 %	> 99.5 %
Compatibility with Gypsum	Pass	Pass
Detail Reproduction	Pass	Pass

## Composition

Virtual impression materials are addition-reaction silicones containing vinylpolysiloxane, methylhydrogensiloxane, organoplatinic complex, silicon dioxide and food dyes.

## Indication

The Virtual line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impression used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from “wax-ups” or for the treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

## Contraindication

The use of Virtual impression materials is contraindicated if the patient is known to be allergic to any of the ingredients in Virtual impression materials.

## Side effects

None known to date.

## Interactions

### The setting of vinylpolysiloxanes is inhibited by latex gloves.

Do not touch preparations or retractions cords with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities, specifically when hand mixing putty. Other materials which may inhibit the set of impression materials include rubber dam, retraction cords and retraction cord fluids. If the operator suspects that the preparation has been contaminated, it is recommended that the preparation be rinsed and dried to eliminate all traces of impurities.

## Application

- **Heavy Body (blue):** High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the tray/wash impression technique.
- **Monophase (blue):** Medium viscosity, hydrophilic vinylpolysiloxane for final impressions. This material is recommended for use as a tray and wash material when using the double-arch, single-phase impression technique.
- **Instructions for Virtual 380 cartridge materials:**  
The Virtual 380 cartridge is suitable for use in conjunction with the Pentamix®, Pentamix® 2 and MixStar mixing units. Carefully remove the Virtual 380 cartridges from the packaging. CAUTION: Avoid dropping the heavy Virtual 380 cartridges, as this may severely damage the cartridges, so that they can no longer be used. Pentamix®, Pentamix® 2 and MixStar are not registered trade marks of Ivoclar Vivadent.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

#### – Inserting the Virtual 380 cartridge:

(Please also refer to the corresponding directions for use of the manufacturer of the mixing unit)

1. Hold the cartridge with the sealed outlet openings facing upwards.
2. Grasp the flap on the sealing cap and bend it upward 90 degrees. Do not apply excessive force in the process.
3. While holding the cartridge securely in one hand, take the flap between thumb and index finger and completely remove the sealing cap. (Fig. 1)
4. Insert the cartridge into the mixing unit according to the manufacturer's instructions. (Fig. 2)
5. To control the flow of the pastes, start the mixing unit according to the manufacturer's instructions and wait until a small amount of material is extruded (This should always be done prior to starting the mixing process). Carefully wipe off excess material in a vertical direction to avoid that the two components become intermixed inside the openings. (Fig. 3)
6. Attach a new mixing tip to the cartridge following the instructions of the mixing unit manufacturer. (Fig. 4)

(NOTE: Only use Ivoclar Vivadent mixing tips – As the mixing tips of other manufacturers do not fit the openings precisely, leakage may occur.)

7. Should you have difficulty in attaching the mixing tip, check whether the central hexagonal socket at the lower end of the tip is properly aligned with the hexagonal drive shaft.
8. Once the mixing tip is correctly seated, place the yellow locking ring over the mixing tip and push it down completely. Turn the ring 1/4 turn clockwise to lock the mixing tip securely in place. (Fig. 5)
9. Load a separate intra-oral elastomer syringe (only Virtual 380 Monophase) or a suitable impression tray with the mixed Virtual impression material.
10. Immediately after extrusion of the material, fully retract the plungers of the mixing unit to avoid continual pressure on the cartridge and prevent the escape of excess material through the mixing tip.
11. Leave the used mixing tip on the cartridge. It serves as a seal until next use. Additionally, cross-contamination of base and catalyst and clogging of the outlet openings with set material is avoided.

**NOTE:** It is recommended that the cartridge be removed from the mixing unit at the end of each working day or when the mixer is not used for a longer period of time. It should be stored in an upright position, with the outlet openings (sealed either by a mixing tip or the sealing cap) facing upwards.

#### – Pre-treatment of impression trays (tray adhesives)

It is highly recommended that tray adhesives (eg: Virtual Tray Adhesive) be used to reduce the chance of distortion when removing impressions from the mouth.

## – Instructions for Virtual Tray Adhesive

1. Make sure that all surfaces of the impression tray are oil-free, clean and dry.
2. Apply a thin layer of Virtual Tray Adhesive to all surfaces of the impression tray (metal or plastic) that will come into contact with the impression material.
3. Allow tray adhesive to dry (approx 3 minutes).
4. Replace bottle cap immediately.

**Note the warnings in the Virtual Tray Adhesive Instructions for Use!**

## – Disinfection of impressions

Impressions made with Virtual impression materials can be immersed in a disinfectant solution (glutaraldehyde 0.5% - benzalkonium chloride 0.5%) immediately for 10 minutes without distortion.

## – Pouring up models

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days. Virtual impression materials are compatible with all popular dental plasters on the market, eg Type 3: Modano® (Heraeus Kulzer), Type 4: Fujirock® (G.C. International).

## – Galvanization

Virtual impression materials can be silver or copper plated in a galvanic bath.

## Special notes

Virtual impression materials should be at room temperature (23 °C / 73.4 °F) when processed. Lower temperatures, eg when stored in the refrigerator, will prolong the time the impression material should remain in the mouth, while higher temperatures will reduce it. Vinylpolysiloxanes are chemically resistant. Unpolymerized materials may stain clothing.

## Warnings

If uncured materials come into contact with the eyes, rinse with copious amounts of water. If irritation persists, seek medical attention. In case of contact with the skin, wash affected areas with plenty of water.

## Storage

- Storage temperature: 2–28 °C / 36–82 °F
- Shelf life: See expiration date on label and packaging
- Keep away from direct heat sources!

**Keep out of the reach of children!**

**For use in dentistry only!**

**Date of preparation:** 08/2011 Rev. 2

**Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein**

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**