Bruksanvisning Brugsanvisning Käyttöohjeet Bruksanvisning Productinformatie Οδηγίες Χρήσεως Kullanım Kılavuzu Инструкция

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Instructions for Use

Description

SpeedCEM is a self-adhesive, self-curing dental resin cement with light-curing option for the cementation of indirect restorations made of high-strength all-ceramic, composite, metal-ceramic and metal. When SpeedCEM is used, additional dentin bonding agents are not required. The automix syringe enables a time-saving direct application of the luting material into the restoration.

Shades

SpeedCEM is available in three shades with various degrees of translucency:

- transparent
- yellow
- white opaque

Working and setting times

The following times apply as soon as SpeedCEM has been dispensed from the automix syringe. (The working and setting times depend on the ambient temperature):

	At room temperature 23 °C ± 1 °C	Intra-orally approx. 37 °C ± 1 °C
Working time	150 ± 30 sec	110 ± 30 sec
Setting time (without working time)	240 ± 60 sec	160 ± 40 sec

Attention

After SpeedCEM has been dispensed from the automix syringe, the setting process can be accelerated by intensive illumination (operating light, ambient light).

Composition

The monomer matrix is composed of dimethacrylates and acyclic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are catalysts, stabilizers and colour pigments (< 1%). The primary particle size of the inorganic fillers is between 0.1 µm and 7 µm. The mean particle size is 5 µm. The total content of inorganic fillers is approx. 40 vol%.

Indication

Permanent cementation of indirect restorations made of

- metal and metal-ceramic (inlays, onlays, crowns, bridges, endodontic posts)
- high-strength all-ceramic; zirconium oxide, lithium disilicate and aluminium oxide ceramics (crowns, bridges, endodontic posts)
- fibre-reinforced composite (endodontic posts)

Contraindication

SpeedCEM is contraindicated

- in situations where the preparation does not provide adequate retention (e.g. veneers).
- if a patient is known to be allergic to any of the ingredients of SpeedCEM.
- in general, SpeedCEM should not be applied to exposed pulp or to dentin which is close to the pulp.

Side effects

Systemic side effects are not known to date. In rare cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol) inhibit polymerization. Consequently, the application of materials containing these substances is to be avoided. Only use eugenol-free temporary cements (e.g. Systemp®.link, Systemp.cem), as eugenol may inhibit the polymerization of SpeedCEM.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Thus the tooth substance must not be disinfected using oxidative agents. Alkaline jet media (e.g. Cojet) impair the effect of SpeedCEM.

Application

1. Removal of the temporary restoration and thorough cleaning of the cavity

Remove possible residue of the temporary luting composite from the cavity or preparation with a polishing brush and an oil- and fluoride-free cleaning paste. Rinse with water spray.

Then lightly dry with water- and oil-free air. Avoid overdrying.

Note:

Cleaning with alcohol can lead to overdrying of the dentin.

2. Try-in of the restoration

Next, the shade, the accuracy of fit and the occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make corrections with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces. It is essential to ensure reliable isolation of the operative area with e.g. OptraGate, cotton rolls or saliva ejector. Dental hard tissue contaminated by saliva needs to be cleaned again as indicated under point 1.

3. Cleaning and pre-treatment of the restoration

3.1 High-strength glass-ceramics (e.g. IPS e.max® Press, IPS e.max CAD)

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 20 seconds (IPS e.max Press, IPS e.max CAD). For all other products, observe the instructions of the manufacturer!
- Thoroughly rinse the restoration with water spray and dry with oil-free air.
- Apply Monobond Plus or any other silanizing agent to the pre-treated surfaces with a brush or microbrush and allow it to react for 60 seconds. Subsequently, disperse Monobond Plus with a strong stream of air.

3.2 Clean and pre-treat the restorations made of other materials according to the instructions of the manufacturer.

Note:

- **IMPORTANT!** In order to achieve a reliable bond, do not clean zirconium oxide and base metal surfaces with phosphoric acid. The surfaces can be cleaned with alcohol.
- Contamination with saliva, blood or water must absolutely be avoided during the entire procedure.
- In case of contamination, the entire pre-treatment of the restoration is to be repeated (see section 3).

4. Application of SpeedCEM into the restoration

For each application, place a new automix tip on the double syringe. Extrude SpeedCEM from the automix syringe and apply the desired quantity directly into the restoration. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until needed once again (replace with a new tip just before the next use).

5. Placement of the restoration and removal of excess cement

a) solely self-curing

Seat the restoration and retain it in place exerting uniform pressure. Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins).

b) self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridges)

Seat the restoration and fix/hold it in place exerting uniform pressure. Excess material is light-cured with the polymerization light (approx. 650 mW/cm², e.g. bluephase®, LOP mode,) for 1 second per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal) at a distance of approx. 0-10mm. Subsequently, it can be easily removed with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins). After complete curing, the removal of excess can be difficult in these areas. Subsequently, light-cure all cement joints again for 20 seconds (approx. 1,200 mW/cm², e.g. bluephase, HIP mode) for final polymerization.

Note

As with all composites, SpeedCEM is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

6. Finishing of the completed restoration

Apply finishing and polishing strips in the proximal regions.

- Check occlusion and functions and adjust if necessary.
- Polish restoration margins with silicone polishers (OptraPol, Astropol®, Politip®-F) or discs.

Special note

SpeedCEM should have room temperature when it is applied.
Cool temperatures can render the material difficult to extrude and increase the processing and curing times.

Warning

Unpolymerized SpeedCEM paste is an irritant. Avoid contact with the skin, mucous membrane and eyes. If SpeedCEM comes into contact with the eyes, immediately rinse with copious amounts of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Do not use SpeedCEM after the indicated date of expiration.
- Cool storage (2–8 °C/ 36–46 °F) for SpeedCEM.
- Store SpeedCEM automix syringe with mixing tip attached after use.
- Shelf life: see expiration date.

Keep material out of the reach of children!

For use in dentistry only!

Date information prepared: 02/2010, Rev 1

Manufacturer:

Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.