Virtual[®] XD

VPS Impression Materials

Instructions for Use Mode d'emploi Instrucciones de Uso

- Vinylpolysiloxane impression material
- Matériau d'empreinte à base de vinyle polysiloxane
- Material de impresión de Vinilpolisiloxano

R_x ONLY

648982 Rev 3 6/2020

ivoclar
vivadent

Virtual[®] XD

VPS Impression Materials

English

Instructions for Use

Description

Virtual XD impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual XD impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

Colors

See table "Technical Data"

Working Time and Setting Times

Material	Speed	Total Work Time [min:sec]	Mouth Set Time* [min:sec]
Extra Light Body [wash]	Regular Set	1:25	3:30
	Fast Set	1:15	2:30
Light Body [wash]	Regular Set	1:25	3:30
	Fast Set	1:15	2:30
Heavy Body [tray]	Regular Set	1:25	4:30
	Fast Set	1:15	2:30
Putty [tray]	Regular Set	1:25	4:30
	Fast Set	1:15	2:30

– Set times are (+/-) 15 seconds

* Minimum time the Impression Material should remain in the mouth before removal.

Technical Data

	Extra Light Body	Light Body	Heavy Body	Putty
Color	Yellow	Orange	Blue	Blue
Classification ISO 4823 / ADA Spec. No. 19	Type 3 light- bodied consistency	Type 3 light- bodied consistency	Type 1 heavy- bodied consistency	Type 0 putty consistency
Mixing Ratio [Base:Catalyst]	1:1	1:1	1:1	1:1

Composition

Virtual XD impression materials are addition-reaction silicones containing vinylpolysiloxane, methylhydrogensiloxane, organoplatinic complex, silica and food dyes.

Indication

The Virtual XD line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impression used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from “wax-ups” or for treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

Contraindication

The use of Virtual XD impression materials is contraindicated if the patient is known to be allergic to any of the ingredients in Virtual XD impression materials.

Side Effects

None known to date

Interactions

The setting of vinylpolysiloxanes is inhibited by latex gloves.

Do not touch preparations or retraction cords with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities, specifically when hand mixing putty. Other materials which may inhibit the set of impression materials include rubber dam, retraction cords and retraction cord fluids. If the operator suspects that the preparation has been contaminated, it is recommended that the preparation be rinsed and dried to eliminate all traces of impurities.

Application

– Putty (Regular and Fast Set)

Very high viscosity vinylpolysiloxane that provides easy mixing and the hydraulic forces necessary to extend wash material into the sulcus providing a detailed impression of prepared margins. The putty viscosity maintains outstanding detail reproduction when used alone. Recommended impression techniques include Putty/Wash One-Step and Putty/Wash Two-Step.

Important: Avoid wearing latex gloves when handling putty, putty jars, and dosing spoons (see note under interactions).

1. Take equal amounts of base (blue) and catalyst (white) using the color coded dosing spoons.
2. Hand mix the identical portions of Virtual XD Putty base and catalyst until you obtain an evenly colored mixture (approximately 30 seconds of mixing time). A proper mixture should have no streaks. Note that equal proportions of base and catalyst must be used in order to obtain proper setting times. A larger quantity of catalyst will not accelerate the setting time.
3. Place mixture in impression tray. It is recommended that the impression trays are pre-coated with Virtual tray adhesive.

– Extra Light Body, Light Body, Heavy Body (Regular and Fast)

Extra Light-Body (Yellow): Extra low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Light Body (Orange): Low viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a wash material when using the Tray/Wash and Putty/Wash impression techniques.

Heavy Body (Blue): High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the Tray/Wash impression technique.

– Instructions for Cartridge Materials

Loading the Cartridge Dispenser

(These instructions are written for use with the Virtual manual dispenser. See separate manufacturers instructions for use if using other dispensers)

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

1. Press the black release lever located below the plunger on the back of the dispenser and pull plunger as far back as possible. (Fig. 1)
2. Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock. (Fig. 2)
3. Once the cartridge is secure in the dispenser, remove the cartridge cap by turning it a 1/4 turn counter clockwise. Discard the cap. (Fig. 3)
4. It is VERY IMPORTANT to clear or bleed the cartridge PRIOR to applying the mixing tip. This will assure a proper mix and set time. Gently press on the dispenser lever until both components (base and catalyst) begin to flow out of the cartridge evenly onto a mixing pad.
5. Select the appropriate mixing tip:
 - a. The large mixing tips (BLUE colored base) are recommended for Heavy Body impression materials
 - b. The small mixing tips (YELLOW colored base) are recommended for Extra-Light Body and Light Body impression materials(The color of the mixing tip base and the color of the cartridge cap are identical)
Insert the mixing tip in a clean motion so that the "V" shaped notch on the colored base is aligned with the "V" shaped notch on the housing of the cartridge. While gripping the colored base, and not the mixing tip, turn the mixing tip base 1/4 turn clockwise. (Fig. 4)
6. If necessary apply an intra-oral tip to the mixing tip (Fig. 5):
 - a. The large intra-oral tips (CLEAR) are to be used with the large mixing tips (BLUE base)
 - b. The small intra-oral tips (YELLOW) are to be used with the small mixing tips (YELLOW base)
7. Begin mixing and dispensing by pressing gently on the lever of the dispenser.

Removing the cartridge

1. Lift the release lever and pull plunger back. Lift the cartridge lock and remove cartridge.
2. Leave the mixing tip on the cartridge. The material will polymerize within the mixing tip, functioning as a natural cap.

– Pre-Treatment of Impression Trays (Tray Adhesives)

It is highly recommended that a tray adhesive (eg: Virtual tray adhesive) be used to reduce the chance of distortion when removing impressions from the mouth.

Instructions for Virtual tray adhesive:

1. Be sure that all surfaces of the impression tray are oil-free, clean and dry.
2. Apply a thin layer of Virtual tray adhesive using the brush provided to all surfaces of the impression tray (metal or plastic) that will come into contact with the impression material.
3. Allow tray adhesive to dry (approximately 3 minutes).
4. Replace bottle cap immediately.

– **Note the warnings in the Virtual tray adhesive Instructions for Use!**

– **Disinfection of Impressions**

Impressions made with Virtual XD impression materials can be immersed in a disinfection solution (glutaraldehyde 0.5% - benzalkonium chloride 0.5%) immediately for 10 minutes without distortion.

– **Pouring Up Models**

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days. Virtual XD impression materials are compatible with all popular dental plasters on the market, e.g. Type 3: Moldano® (Heraeus Kulzer), Type 4: Fujirock® (G.C. International)

– **Galvanization**

Virtual XD impression can be silver or copper plated in a galvanic bath.

Special Notes

Vinylpolysiloxanes are chemically resistant. Unpolymerized materials may stain clothing.

Warnings

If uncured materials come into contact with the eyes, rinse with copious amounts of water. If irritation persists seek medical attention. In case of contact with the skin, wash affected areas with soap and water.

Storage

- Storage temperatures: 2–28 °C / 36–82 °F
- Shelf life: See expiration date on label and packaging
- Keep away from direct heat sources!

Keep out of the reach of children! For use in dentistry only!

Date information prepared: 3/2013 Rev.2

Ivoclar Vivadent, Inc. Amherst, New York, USA

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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©2020 Ivoclar Vivadent, Inc. All rights reserved.
trademarks of Ivoclar Vivadent.
Ivoclar Vivadent Inc. Amherst, New York, USA