ivoclar
vivent

Short and long span bridges, PFM crowns, posts. Recommended for inlays, onlays, 3/4 crowns, telescopes, crowns, conus crowns, posts.

Composition

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Yellow, gold-based dental metal ceramic casting alloy, type 3

ivoclar Vivadent Worldwide

Brazil	India	Russia
ivoclar Vivadent Ltda. Rua Geraldo Flausino Gomes, 78 - 6.º andar Cxj. 61162 Bairro: Brooklin Novo Veria Desai Road, Andheri (West) Mumbai, 400 053 Tel. +91 (22) 2673 0302 Fax +91 (22) 2673 0301	ivoclar Vivadent Marketing Ltd 503/504 Raheja Plaza 15 B Shah Industrial Estate Veer Desai Road, Andheri (West) Mumbai, 400 053 Tel. +91 (22) 2673 0302 Fax +91 (22) 2673 0301	ivoclar Vivadent Marketing Ltd. Derbenenskaja Naberenskaja 11, Geb. V 115114 Moscow Tel. +7 495 913 66 19 Fax +7 495 916 66 15
USA	Canada	Singapore
ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.us	ivoclar Vivadent AG Benderstrasse 2 FL-9494 Schaan Tel. +423 235 35 35 Fax +423 235 33 60 www.ivoclarvivadent.com	ivoclar Vivadent Marketing Ltd. 171 Chin Swee Road Singapore 169877 Tel. +65 6535 6175 Fax +65 6535 4991
Australia	China	Japan
ivoclar Vivadent Pty. Ltd. 1 - 5 Overseas Drive P.O. Box 367 Noble Park, Vic. 3174 Tel. +61 3 979 595 99 Fax +61 3 979 596 45 www.ivoclarvivadent.com.au	ivoclar Vivadent Marketing Ltd. Rm 602 Kuen Tang International Business Plaza No. 798 Zhao Jia Bang Road Shanghai 200030 Tel. +86 21 5456 0776 Fax +86 21 6445 1561 www.ivoclarvivadent.com.cn	ivoclar Vivadent K.K. 1-28-24-4F Horago Bunkyo-ku, Tokyo 113-0033 Tel. +81 3 6903 3535 Fax +81 3 5844 3657 www.ivoclarvivadent.co.jp
Austria	Mexico	Spain
ivoclar Vivadent GmbH Bremschtr. 16 F-74479 Saint-Jorioz Tel. +43 5552 624 49 Fax +43 5552 624 49 www.ivoclarvivadent.com	ivoclar Vivadent S.A. de C.V. Av. Mazatlan No. 61, Piso 2 Col. Condesa 06170 México City, D.F. Bogotá Tel. +52 (55) 50 62 10 29 Tel. +52 (55) 50 62 10 29 www.ivoclarvivadent.com.mx	ivoclar Vivadent S.A. c/ Emilio Muñoz Nº 15 E-28037 Madrid Tel. +34 913 75 78 20 Fax +34 913 75 78 38 www.ivoclarvivadent.com
UK	Turkey	Sweden
ivoclar Vivadent UK Limited Ground Floor Compass Building Fieldspar Close Warrens Business Park Enderby Leicester LE19 4SE England Tel. +44 (0) 79 61 17 88-0 Fax +44 (0) 79 61 17 88-0 www.ivoclarvivadent.co.uk	ivoclar Vivadent Liason Office Ahi Evran Caddesi no 1 Polaris Is Merkezi Kat: 7 80670 Maslak Istanbul Tel. +90 536 577 12 62	ivoclar Vivadent AB Dahlgrän 14 S-169 56 Solna Tel. +46 (0) 8 514 93 930 Fax +46 (0) 8 514 93 940
Germany	Poland	USA (Continental)
ivoclar Vivadent GmbH Dr. Adolf-Schneider-Str. 2 D-73479 Ellwangen, Jost Tel. +49 (0) 79 61 17 88-0 Fax +49 (0) 79 61 17 88-0 www.ivoclarvivadent.de	ivoclar Vivadent Sp. z o.o. ul. Jana Pawła I 78 PL-00175 Warszawa Tel. +48 22 635 54 96 Fax +48 22 635 54 69 www.ivoclarvivadent.pl	ivoclar Vivadent Inc. 175 Pineview Drive Amherst, NY 14228 Tel. +1 800 533 6825 Fax +1 716 691 2285 www.ivoclarvivadent.us

Instructions for Use

Instructions de use

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Brugsanvisning

Modo d'empio

Modélado en cera

Modélado en cera

Modélado en cera

ISO 13485
Quality Management System Certified

MADE IN U.S.A.
unless otherwise indicated
Made in U.S.A.
(fabrique au Canada)

ivoclar vivadent
passion vision innovation

Rev. Date 1/10
Rev. 3

Caution: US Federal Law restricts this device to sale by or on the order of a licensed dentist.

LOT

EN INSTRUCTIONS FOR USE

MODELLATION
Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

BURN-OUT
The suggested burnout temperature: 750-820C/1380-1510F

MELTING AND CASTING
Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

Casting Temperature: 1260-1320C/2300-2410F

METAL PREPARATION
Carefully divest and clean the object with AL₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

OXIDATION
Blast surface with 50-100 micron AL₂O₃ at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 1010C/1850F without vacuum with 5 min. hold time at top temperature.

HEAT TREATMENT
Hardening: 540C/1000F for 10 minutes; air cool.

SOLDERS AND FLUXES
Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Pre-Solder: High Fusing Yellow Ceramic Solder, High Fusing Bondal Flux
Post-Solder: .615 Fine Solder, Low Fusing White Gold Solder, .585 Fine Solder, Bondal Flux
Laser weld material: Laser Ceramic Yellow

POLISHING
After glazing remove oxide and flux residue and finish/polish the framework with rubber finish-ers/polishers.

INDICATIONS
Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

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