**ivoclar vivadent**

ISO 13485  
  

Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

**Indications**

Indicated for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

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**Yellow, gold-based dental metal ceramic casting alloy, Type 4**

**Produktformate**  
 Instruktions for Use  
 Gebrauchsinformation  
 Instruções de Uso  
 Instrucciones de Uso  
 Mode d'emploi  
 Brugsanvisning  
 Instruções de Uso  
 Brugsanvisning  
 Instruções de Uso  
 Käyttöohjeet

**Brugsanvisning**

**Mode d'emploi**

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**Mode d'emploi**

**Brugsanvisning**

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**ivoclar Vivadent Worldwide**

**Brazil**  
 Ivoclar Vivadent Ltda.  
 Rua Geraldo Flausino Gomes,  
 15 B Shah Industrial Estate,  
 6.º andar, Cj. 61/62  
 Bairro: Brooklin Novo  
 CEP: 04575-060 - São Paulo - SP  
 Tel: +55 (11) 3466 0800  
 Fax: +55 (11) 3466 0840  
 www.ivoclarvivadent.com.br

**Canada**  
 Ivoclar Vivadent Inc.  
 2785 Skyway Avenue, Unit 1  
 Mississauga, Ontario, L4W4Y3  
 Tel: +1 905 238 5700  
 Fax: +1 905 238 5711  
 www.ivoclarvivadent.ca

**China**  
 Ivoclar Vivadent Marketing Ltd.  
 Rm 603 Xuen Tang  
 International Business Plaza  
 No. 798 Zhao Jia Bang Road  
 Shanghai 200030  
 Tel: +86 21 5456 0776  
 Fax: +86 21 6445 1561  
 www.ivoclarvivadent.com

**Colombia**  
 Ivoclar Vivadent Marketing Ltd.  
 Calle 134 No. 7-83, Of. 520  
 Bogotá  
 Tel: +57 (1) 627 33 99  
 Fax: +57 (1) 633 16 63  
 www.ivoclarvivadent.com

**France**  
 Ivoclar Vivadent SAS  
 B.P. 118  
 F-74413 Saint-Jorioz  
 Tel: +33 450 88 64 00  
 Fax: +33 450 88 61 52  
 www.ivoclarvivadent.fr

**Austria**  
 Ivoclar Vivadent GmbH  
 Bremschtr. 16  
 Postfach 223  
 A-6706 Burs  
 Tel: +43 5552 624 49  
 www.ivoclarvivadent.com

**Australia**  
 Ivoclar Vivadent Pty. Ltd.  
 1 - 5 Overseas Drive  
 P.O. Box 367  
 Noble Park, Vic. 3174  
 Tel: +61 3 979 595 99  
 Fax: +61 3 979 596 45  
 www.ivoclarvivadent.com.au

**Germany**  
 Ivoclar Vivadent GmbH  
 Dr. Adolf-Schneider-Str. 2  
 D-73479 Ellwangen, Germany  
 Tel: +49 (0) 79 61 1 889-0  
 Fax: +49 (0) 79 61 63 26  
 www.ivoclarvivadent.de

**USA**  
 Ivoclar Vivadent Inc.  
 175 Pineview Drive  
 Amherst, NY 14228  
 Tel: +1 800 533 6825  
 Fax: +1 716 691 2285  
 www.ivoclarvivadent.us

**India**  
 Ivoclar Vivadent Marketing Ltd.  
 503/504 Raheja Plaza  
 Veera Desai Road, Andheri (West)  
 Mumbai, 400 053  
 Tel: +91 (22) 2673 0302  
 Fax: +91 (22) 2673 0301  
 www.ivoclarvivadent.com

**Italy**  
 Ivoclar Vivadent s.r.l. & C. s.a.s.  
 Via Gustav Flora 32  
 39025 Naturno (BZ)  
 Tel: +39 0473 67 01 11  
 Fax: +39 0473 66 77 80  
 www.ivoclarvivadent.it

**Spain**  
 Ivoclar Vivadent S.A.  
 c/ Emilio Muñoz Nº 15  
 E-28037 Madrid  
 Tel: +34 913 75 78 20  
 Fax: +34 913 75 78 38  
 www.ivoclarvivadent.com

**Sweden**  
 Ivoclar Vivadent AB  
 Dalvägen 14  
 S-169 56 Solna  
 Tel: +46 (0) 8 514 93 930  
 Fax: +46 (0) 8 514 93 940  
 www.ivoclarvivadent.se

**Turkey**  
 Ivoclar Vivadent Liaison Office  
 Ahi Evran Caddesi no 1  
 Polaris is Merkezi Kat: 7  
 08070 Maslak  
 Istanbul  
 Tel: +90 536 577 12 62

**UK**  
 Ivoclar Vivadent UK Limited  
 Ground Floor Compass Building  
 Feldspar Close  
 Warrens Business Park  
 Enderby  
 Leicestershire LE19 4SE England  
 Tel: +44 116 284 78 80  
 Fax: +44 116 284 78 81  
 www.ivoclarvivadent.co.uk

**Manufacturer**

**USA**  
 Ivoclar Vivadent Inc.  
 175 Pineview Drive  
 Amherst, NY 14228  
 Tel: +1 800 533 6825  
 Fax: +1 716 691 2285  
 www.ivoclarvivadent.us

**ISO 13485**  
 Quality Management System Certified

**Made in U.S.A. unless otherwise indicated**  
 Made in Canada (fabrique au Canada)

**www.ivoclarvivadent.com**

**Rev. Date 1/10**  
 Rev. 7

**Caution:** US Federal Law restricts this device to sale by or on the order of a licensed dentist.

**LOT**

**ivoclar vivadent**  
 passion vision innovation

**EN INSTRUCTIONS FOR USE**

**MODELLATION**  
 Design the framework in a reduced anatomical shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

**SPRUNG**  
 Provide the modelled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

**INVESTMENT**  
 Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

**BURN-OUT**  
 The suggested burnout temperature: 750-820C/1380-1510F

**MELTING AND CASTING**  
 Use a separate carbon/ceramic crucible for each alloy and preheat the crucible (ceramic) in the burnout furnace. Used and new alloys must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. If you use Ivoclar Vivadent Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting bench cool to room temperature.

**Casting Temperature:** 1260-1320C/2300-2410F

**METAL PREPARATION**  
 Carefully divest and clean the object with Al<sub>2</sub>O<sub>3</sub>. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

**OXIDATION**  
 Blast surface with 50-100 micron Al<sub>2</sub>O<sub>3</sub> at max. 4.5 bar/65 psi pressure before oxidation. Subsequently, clean with ultrasonic or steam. Place the object on the firing tray and provide adequate support. Place the tray in a porcelain furnace set at a low temperature of 650C/1200F and increase the furnace temperature to 1010C/1850F without vacuum with 5 min. hold time at top temperature.

**HEAT TREATMENT**  
 Hardening: 510C/950F for 30 minutes; air cool.

**SOLDERS AND FLUXES**  
 Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

**Pre-Solder:** High Fusing Yellow Ceramic Solder, High Fusing Bondal Flux  
**Post-Solder:** .615 Fine Solder, Low Fusing White Gold Solder, .585 Fine Solder, Bondal Flux  
**Laser weld material:** Laser Ceramic Yellow

**POLISHING**  
 After glazing remove oxide and flux residue and finish/polish the framework with rubber finish-ers/polishers.

**INDICATIONS**  
 Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, posts, short and long span bridges, PFM crowns.

**CONTRAINDICATIONS**  
 For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consult with a physician is recommended.

**SIDE EFFECTS**  
 In individual cases, sensitivity or allergies to components of this alloy may occur.

**INTERACTIONS**  
 Galvanic effects may occur between different alloys in the same oral environment.  
 For additional information look into the alloy property chart.

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