ivoclar vivadent

Indications

crowns, short and long bridges, posts.
 crowns, 3/4 crowns, telescope crowns, conus crowns, short and long bridges, posts.

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Yellow, gold-based dental casting alloy, Type 3

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INSTRUCTIONS FOR USE

MODELLATION
 Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration. Single crowns require a thickness of minimum 0.3 mm; abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used. If a composite or resin veneer is required, mechanical retention is recommended.

SPRUNG
 Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center. The connection sprues between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
 Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instruction.

BURN-OUT
 The suggested burnout temperature:
 High heat temperature investment: 650-760C/1200-1400F
 Low heat temperature investment: 480-540C/900-1000F

MELTING AND CASTING
 Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturers instructions for use. Ideally a compressed air and natural gas torch should be used to melt C&B alloys because propane and oxygen is much too hot and can easily overheat these alloys. If you are using propane and oxygen the pressure should be a 0.15 bar/2 psi for propane and 0.35 bar/5 psi for oxygen. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Use casting flux if needed. After casting bench cool to room temperature.

METAL PREPARATION
 Carefully divest and clean the object with AL₂O₃ glass beads, or a pickling agent (such as Prevoxx®). Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Avoid inhalation of dust during grinding!

HEAT TREATMENT
 Annealing: 620C/1150F for 10 minutes; quench immediately (water)
 Hardening: 345C/650F for 10 minutes; air cool.

SOLDERS AND FLUXES
 Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

Polishing
 After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS
 Recommended for inlays, onlays, 3/4 crowns, crowns, telescope crowns, conus crowns, short and long span bridges, posts.

CONTRAINDICATIONS
 For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
 In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
 Galvanic effects may occur between different alloys in the same oral environment.

For additional information look into the alloy property chart.

SPRUNG METHOD

RECOMMENDATIONS

DIRECT: single crowns, inlays and onlays

INDIRECT: multiple units and multiple single crowns

INSTRUCTIONS:

- Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.
- Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.
- Make sure that the sprue is connected to the thickest part of the restoration.
- The connection between the sprue and the restoration should be flared ("trumpet" configuration); eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification.
- Maintain proper alloy weight in order to eliminate the negative effect of a too large button during solidification. The rule of thumb for proper alloy weight determination is: total wax weight x specific gravity of the alloy = proper alloy weight.

SUGGESTIONS:

- Thickness and design of wax pattern: follow your manufacturers instructions.
- Use chill-vents when casting heavy and/or large restorations.

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