Indications: Inlays, onlays, 3/4 crowns, crowns, telescopic crowns, conus crowns, short and long bridges, posts.

Table with 10 columns: Alloy, Weight, Density, etc. for various dental alloys.

Composition


Instructions for Use: Production of wax patterns, casting, and finishing instructions for various dental applications.

White, silver-based dental casting alloy, Type 4

Worldwide distribution list with columns for Country, Address, Phone, and Website for various regions including Brazil, India, Russia, Canada, Italy, China, USA, Colombia, Australia, Austria, Germany, and Japan.

ISO 13485 Quality Management System Certified. Made in U.S.A. unless otherwise indicated. Made in Canada (fabrique au Canada).

EN INSTRUCTIONS FOR USE

MODELLATION Wax to full contour for crown and bridge frame design. For composite, build up the framework in a reduced anatomic shape taking the planned veneer into consideration.

SPRUNG Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or indirect technique being sure that the reservoir is positioned in the heat center.

INVESTMENT Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used.

BURN-OUT The suggested burnout temperature: 500-700C/940-1300F. High heat temperature investment: 480-540C/900-1000F.

MELTING AND CASTING Use a separate carbon/ceramic crucible for each alloy. Used and new alloy must be in a ratio of 1:1. Depending on the type of casting machine, follow the manufacturer's instructions for use.

METAL PREPARATION Carefully divest and clean the object with Al2O3 glass beads, or a pickling agent (such as Prevoxx®).

HEAT TREATMENT Annealing: 705C/1300F for 15 minutes; quench immediately (water). Hardening: 400C/750F for 30 minutes; air cool.

SOLDERS AND FLUXES Design the soldering patty as small as possible and preheat it in the furnace at approximately 600C/1112F.

POLISHING After soldering or heat treatment, remove oxide and flux residue and finish and polish the framework with rubber finishers and polishers.

INDICATIONS Recommended for inlays, onlays, 3/4 crowns, crowns, telescopic crowns, conus crowns, short and long span bridges, posts. CONTRAINDICATIONS For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SPRUNG METHOD

RECOMMENDATIONS DIRECT: single crowns, inlays and onlays. INDIRECT: multiple units and multiple single crowns.

INSTRUCTIONS: 1. Select a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration. 2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.

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