# Virtual<sup>®</sup> XD 380

VPS Impression Materials

**Instructions for Use**

**Mode d'emploi**

**Istruzioni d'uso**

- Vinylpolysiloxane impression material
- Matériau d'empreinte à base de vinyle polysiloxane
- Materiale per impronte al polivinilsilossano (silicone d'addizione)

**R<sub>x</sub> ONLY**

655759 Rev. 3 6/2020

**ivoclar**  
**vivadent**  
passion vision innovation

## Instructions for Use

### Description

Virtual XD impression materials are addition-reaction silicones (vinylpolysiloxanes) used to create fine detail impressions of dentition. Virtual XD impression materials are available in a variety of viscosities allowing dental professionals to choose the material and technique best suited for each individual case.

### Colors

See table "Technical Data"

### Working times and setting times

| Material             | Speed       | Total Work Time<br>[min:sec] | Mouth Set Time*<br>[min:sec] |
|----------------------|-------------|------------------------------|------------------------------|
| Heavy Body<br>[tray] | Regular Set | 1:25                         | 4:30                         |
|                      | Fast Set    | 1:15                         | 2:30                         |

– Set times are (+/-) 15 seconds

\* Minimum time the Impression Material should remain in the mouth before removal.

### Technical Data

|   | Heavy Body                         |
|---|------------------------------------|
| Color   | Blue                               |
| Classification<br>ISO 4823 / ADA Spec. No. 19 | Type 1<br>Heavy-bodied consistency |
| Mixing Ratio<br>[Base : Catalyst]             | 5:1                                |

## Composition

Virtual XD impression materials are addition-reaction silicones containing vinylpolysiloxane, methylhydrogensiloxane, organoplatinic complex, silicon dioxide and food dyes.

## Indication

The Virtual XD line of addition-silicone (vinylpolysiloxane) impression materials is recommended for use to create highly detailed impressions of the hard and soft tissues of the oral cavity.

- Final impression used for the fabrication of indirect restorations (crowns, bridges, inlays, onlays and veneers)
- Dental implant impressions
- Matrix from “wax-ups” or for the treatment planning, study models
- Edentulous impressions
- Matrix used to create temporary restorations

Each viscosity provides specific attributes to meet the needs of the dental professional.

## Contraindication

The use of Virtual XD impression materials is contraindicated if the patient is known to be allergic to any of the ingredients in Virtual XD impression materials.

## Side effects

None known to date.

## Interactions

**The setting of vinylpolysiloxanes is inhibited by latex gloves.**

Do not touch preparations or retractions cords with latex gloves. It is recommended that operators wash their hands thoroughly or use vinyl gloves in order to eliminate all traces of impurities, specifically when hand mixing putty. Other materials which may inhibit the set of impression materials include rubber dam, retraction cords and retraction cord fluids. If the operator suspects that the preparation has been contaminated, it is recommended that the preparation be rinsed and dried to eliminate all traces of impurities.

## Application

- **Heavy Body (blue):** High viscosity, hydrophilic vinylpolysiloxane suitable for final impressions. This material is recommended for use as a tray material when using the tray/wash impression technique.
- **Instructions for Virtual XD 380 cartridge materials:**  
The Virtual XD 380 cartridge is suitable for use in conjunction with the dynamic mixing devices. Carefully remove the Virtual XD 380 cartridges from the packaging.  
**CAUTION:** Avoid dropping the heavy Virtual XD 380 cartridges, as this may severely damage the cartridges, so that they can no longer be used.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Fig. 6

#### – Inserting and using the Virtual XD 380 cartridge:

(Please also refer to the corresponding directions for use of the manufacturer of the mixing unit)

1. Hold the cartridge with the sealed outlet openings facing upwards.
2. Grasp the flap on the sealing cap and bend it upward 90 degrees. Do not apply excessive force in the process.
3. While holding the cartridge securely in one hand, take the flap between thumb and index finger and completely remove the sealing cap. (Fig. 1)
4. Insert the cartridge into the mixing unit according to the manufacturer's instructions. (Fig. 2)
5. To control the flow of the pastes, start the mixing unit according to the manufacturer's instructions and wait until a small amount of material is extruded (This should always be done prior to starting the mixing process). Carefully wipe off excess material in a vertical direction to avoid that the two components become intermixed inside the openings. (Fig. 3)
6. Attach a new mixing tip to the cartridge following the instructions of the mixing unit manufacturer. (Fig. 4)  
(NOTE: Only use Ivoclar Vivadent mixing tips – As the mixing tips of other manufacturers do not fit the openings precisely, leakage may occur.)
7. Should you have difficulty in attaching the mixing tip, check whether the central hexagonal socket at the lower end of the tip is properly aligned with the hexagonal drive shaft.
8. Once the mixing tip is correctly seated, place the yellow locking ring over the mixing tip and push it down completely. Turn the ring 1/4 turn clockwise to lock the mixing tip securely in place. (Fig. 5)
9. Load a separate intra-oral elastomer syringe or a suitable impression tray with the mixed Virtual XD impression material.
10. For an optimal mix with the Virtual Mixing Machine select the mixing program marked with the diamond . (Fig. 6).
11. Immediately after extrusion of the material, fully retract the plungers of the mixing unit to avoid continual pressure on the cartridge and prevent the escape of excess material through the mixing tip.
12. Leave the used mixing tip on the cartridge. It serves as a seal until next use. Additionally, cross-contamination of base and catalyst and clogging of the outlet openings with set material is avoided.

**NOTE:** It is recommended that the cartridge be removed from the mixing unit at the end of each working day or when the mixer is not used for a longer period of time. It should be stored in an upright position, with the outlet openings (sealed either by a mixing tip or the sealing cap) facing upwards.

#### – Pre-treatment of impression trays (tray adhesives)

It is highly recommended that tray adhesives (eg: Virtual Tray Adhesive) be used to reduce the chance of distortion when removing impressions from the mouth.

### Instructions for Virtual tray adhesive:

1. Be sure that all surfaces of the impression tray are oil-free, clean and dry.
2. Apply a thin layer of Virtual tray adhesive using the brush provided to all surfaces of the impression tray (metal or plastic) that will come into contact with the impression material.
3. Allow tray adhesive to dry (approximately 3 minutes).
4. Replace bottle cap immediately.

– **Note the warnings in the Virtual tray adhesive Instructions for Use!**

#### – **Disinfection of Impressions**

Impressions made with Virtual XD impression materials can be immersed in a disinfection solution (glutaraldehyde 0.5% - benzalkonium chloride 0.5%) immediately for 10 minutes without distortion.

#### – **Pouring Up Models**

The impression may be poured immediately after disinfection, or up to two weeks later, provided that the impression is stored at room temperature. Dimensional stability is guaranteed for 14 days. Virtual XD impression materials are compatible with all popular dental plasters on the market, e.g. Type 3: Moldano® (Heraeus Kulzer), Type 4: Fujirock® (G.C. International)

#### – **Galvanization**

Virtual XD impression can be silver or copper plated in a galvanic bath.

### **Special Notes**

Vinylpolysiloxanes are chemically resistant. Unpolymerized materials may stain clothing.

### **Warnings**

If uncured materials come into contact with the eyes, rinse with copious amounts of water. If irritation persists seek medical attention. In case of contact with the skin, wash affected areas with soap and water.

### **Storage**

- Storage temperatures: 2–28 °C / 36–82 °F
- Shelf life: See expiration date on label and packaging
- Keep away from direct heat sources!

**Keep out of the reach of children! For use in dentistry only!**

**Date information prepared:** 3/2013 Rev.1

**Ivoclar Vivadent, Inc.** Amherst, New York, USA

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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**[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)**

©2020 Ivoclar Vivadent, Inc. All rights reserved.  
Virtual and Ivoclar Vivadent are registered  
trademarks of Ivoclar Vivadent.  
Rev 3 6/2020

The logo for Ivoclar Vivadent features a series of seven dots of varying sizes arranged in a slight arc above the company name. The word "ivoclar" is in a lowercase, sans-serif font, and "vivadent" is in a larger, bold, lowercase, sans-serif font. A registered trademark symbol (®) is located to the right of "vivadent". Below the company name, the tagline "passion vision innovation" is written in a smaller, lowercase, sans-serif font.

**ivoclar** ·  
**vivadent**®  
passion vision innovation

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**Ivoclar Vivadent Inc.** Amherst, New York, USA