### Bruksanvisning

### Bruksanvisning

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### Bruksanvisning

### Productinformatie

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## Instructions for Use

### Description:

Fermit and Fermit-N are light-curing single-component materials for temporary restorations, generally not requiring the application of a temporary cement.

### Composition:

#### Fermit

The monomer matrix consists of high-molecular dimethacrylates. The fillers are highly dispersed silicon dioxide and copolymers. Catalysts and stabilizers are additional ingredients.

#### Fermit-N

The monomer matrix consists of high-molecular dimethacrylates. The fillers are highly dispersed silicon dioxide and copolymers. Catalysts and stabilizers are additional ingredients.

### Indications:

- Fermit is especially suitable for deep inlay preparations with parallel walls. Fermit can also be used if small undercuts are present.
- Fermit can also be used as a relining material for prefabricated temporary crowns and bridges made of polycarbonate or methacrylate.
- Fermit may also be used for the temporary sealing of implant screw access openings.
- Fermit-N is especially suitable for large preparations (onlays). If undercuts are present, please note that after curing, Fermit-N is less elastic than Fermit, which may complicate removal.

### Contraindications:

Unsuitable as temporary crown and bridge material. Fermit and Fermit-N should not remain in the oral cavity for more than one month. In rare cases, the ingredients in Fermit or Fermit-N may cause a sensitizing reaction. In these cases, the material should no longer be used.

5555840406/ME3/H

**Side effects:**

After longer periods of service, the material may demonstrate discolouration. This, however, does not compromise the performance of the material. The product should not be used in patients with a known allergy to any of its ingredients.

**Interactions:**

Phenolic substances (e.g. eugenol) may inhibit the polymerization of Fermit / Fermit-N. The application of zinc oxide-eugenol cements in combination with Fermit / Fermit-N is thus to be avoided.

**Application procedure:****Preparation appointment**

1. Prepare; place base or liner; take impression as usual.
2. Dry prepared teeth; use a matrix if necessary.
3. Clean cavity with water spray.
4. Application of Systemp.desensitizer:  
Apply Systemp.desensitizer and carefully brush the material into the dentin for 10 seconds.
5. Carefully dry with blown air.
6. Place a sufficient quantity of Fermit or Fermit-N into the cavity with a spatula or another suitable instrument; contour leaving minimum excess.
7. Increments of up to 4 mm may be cured with a standard curing light (e.g. bluephase® C5) in 10 seconds if the light intensity is higher than 500 mW/cm<sup>2</sup>. The same depth of cure is achieved with LED lights that produce a higher light intensity (e.g. bluephase®) after an exposure time of 10 seconds in the low power mode (LOP).
8. Remove matrix
9. Verify occlusion.

**Insertion appointment**

1. Insert a suitable instrument (e.g. probe / scaler) into the temporary restoration. Remove Fermit or Fermit-N from the cavity.
2. Following this, clean the cavity (with e.g. rubber cup and cleaning paste).
3. Try in and seat the restoration.

**Note 1:**

Fermit/Fermit-N generally remains in place in standard inlay preparations. If mechanical retention is minimal Fermit / Fermit-N is cemented with a eugenol-free cement (e.g. Systemp®.link). If Fermit / Fermit-N is applied to undercuts of the proximal region, retention can be improved. In this case, do not place a matrix, but a wooden interdental wedge before applying Fermit / Fermit-N.

**Note 2:**

Because of the similar chemical compositions, Fermit / Fermit-N bonds with light-curing base/liner material. If such materials are used (light-curing glass ionomer cements, composites), isolate the base with glycerine gel (Liquid Strip). This prevents the base/lining from being removed together with Fermit / Fermit-N.

**Note 3:**

Use silicone rubber finishers (Politip-F) or tungsten-carbide finishers for grinding and excess removal. Then polish with silicone rubber polishers (Politip-P). Excess material may also be removed with a scalpel.

**Note 4:**

To increase stability in larger cavities, a composite material (e.g. Heliomolar® or Tetric EvoCeram®) can be applied to the occlusal part of the temporary. Light-cure Fermit / Fermit-N separately or together with the restorative material.

**Note 5:**

Fermit / Fermit-N can be contoured more easily by wetting the instrument with Systemp.desensitizer or an unfilled bonding agent (e.g. Heliobond).

**Warning:**

Unpolymerized Fermit / Fermit-N should not come in contact with skin, mucous membrane, or eyes. Unpolymerized Fermit / Fermit-N may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

**Storage:**

- Do not use Fermit / Fermit-N after the indicated date of expiration.
- Storage temperature (2–28 °C / 36–82 °F)
- Close syringe immediately after use.
- Store material protected from light.
- Shelf life: see date of expiration

**Keep out of the reach of children!****For use in dentistry only!****Date information prepared: 02/2006****Rev. 1****Manufacturer:**

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

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## **Ivoclar Vivadent AB**
## **Ivoclar Vivadent UK Limited**
## **Ivoclar Vivadent, Inc.**