Instructions for Use Gebrauchsinformation

Mode d'emploi

Istruzioni d'uso

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Bruksanvisning

Bruksanvisning

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Bruksanvisning

Productinformatie

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For dental use only.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist.

Made in Liechtenstein

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

Instructions for Use

Description

OptraGate is an auxiliary aid for clinical use which allows the operating field to be accessed with ease. It provides retraction of lips and cheeks during dental treatment and offers the clinician increased visibility and accessibility. The device, which features high flexibility and elasticity in all directions of movement, ensures patient comfort and assists in keeping the patient's mouth open. OptraGate is latex-free.

OptraGate is available in three sizes: Regular, Small and Junior for children between 5 and 10 years of age.

Indication

Suitable for use in diagnostic, preventive and therapeutic dental procedures (in particular restorative procedures, cavity preparation, prophylaxis, diagnostic procedures, scaling & root planing, tooth whitening, impression making).

Contraindications

The use of OptraGate is contraindicated in surgical and implantological procedures.

Description of the procedure

When OptraGate is used, the small inner ring of the device is positioned in the area of the gingivobuccal fold, while the outer ring remains outside the mouth. The elastic component between the two rings thus embraces the lips of the patient and provides retraction due to the restoring force of the rings (see Figure 1).

Fig 1: OptraGate in place. It does not interfere with the patient's occlusion.

1. Select the appropriate sized OptraGate. The size (R for regular, S for small or J for Junior) is printed on the left tab.
2. It is advisable to moisten the OptraGate with water prior to using it, as this will facilitate insertion.
3. The OptraGate may be inserted best from both the 12 o'clock position with the patient lying down, or from the 8 o'clock position with the patient sitting upright (or, the 4 o'clock position if the clinician is left-handed).
4. To ease placement, with the tabs pointing downward, the operator should grasp the inner ring by reaching through the outer ring. The inner ring should be held between the thumb and index finger and pressed together slightly (see Figure 2).

Fig 2: Correct hold on the inner ring with thumb and index finger. The ring is slightly compressed in the process.

5. The slightly compressed inner ring is then inserted into the buccal corridor on one side, so that the elastic component embraces the corner of the mouth and the outer ring is still located extraorally (see Figure 3).

Fig 3: The inner ring is inserted between the teeth and the corner of the mouth.

6. As soon as the OptraGate is securely in place on one side, the other side is inserted in the same manner by slightly bending it (see Figure 4).

Fig 4: The inner ring is positioned behind the two corners of the mouth.

7. Subsequently, the inner ring is placed behind the upper and lower lips and thus OptraGate assumes its final, stable position (see Figure 5).

Fig 5: Final positioning behind the upper and lower lip.

- The securely positioned OptraGate improves the overall view and accessibility of the working field. Lateral movements of the mandible and opening of the mouth to various degrees may provide additional space and access to the treatment area. The OptraGate may even be left in place when checking occlusion.
- For the removal of OptraGate, the outer rim is held and compressed slightly, so that the inner ring is loosened from the gingivobuccal fold. The inner ring is carefully removed from one corner of the mouth. Due to its restoring force it will then slip free on its own.

If, in individual cases, the inner ring should show a tendency to slip out of the mouth in the area of the lower lip when the mouth is closed completely, positioning of the inner ring deeper inside the vestibule is usually sufficient to solve the problem. Sometimes, the use of another size may also be a solution.

Note

OptraGate is for single patient use only.

**Keep out of the reach of children!
For use in dentistry only!**

**Date of preparation: 05/2005
Revision 1**

Manufacturer:

Ivoclar Vivadent AG, 9494 Schaan,
Principality of Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent AB
Ivoclar Vivadent UK Limited
Ivoclar Vivadent, Inc.