

e.max[®] Press
Multi
A2



IPS e.max[®] PRESS MULTI

**THE FIRST-EVER POLYCHROMATIC PRESS INGOT:
PRESS IN MULTICOLOUR – GLAZE – FINISH!**

- Lifelike transition of shade and translucency from dentin to incisal areas
- For the fabrication of monolithic crowns, veneers and hybrid abutment crowns
- Proven IPS e.max Press (LS₂) material for high strength (400 MPa) and longevity
- Coordinated with high-precision Programat press furnaces featuring patented processing technology for utmost efficiency
- Maximum cost effectiveness in the press technique

IPS e.max Press Multi is available in the following shades:
A1, A2, A3, A3.5, B1, B2, C1, C2, D3 and BL2

673701/2014-08/en

www.ivoclarvivadent.com



ivoclar
vivadent[®]
passion vision innovation