

PrograBase X10

The multifunctional software for your Programat® furnace



Operate | Manage | Protocol | Download

PrograBase X10

The multifunctional software for your Programat[®] furnace

The newest version of the Programat computer software now offers even more possibilities.



Operate

Operate the Programat furnace directly from your computer. The screen shows all the important information, such as the Programat number, the firing status and the current firing temperature. This option allows the remote operation and testing of the Programat furnace from your computer/laptop.



Manage

Manage your programs conveniently with your computer. Subsequently, save your acquired data on the furnace and/or your computer. You can also adjust the firing and sintering programs according to your needs, create new programs and copy existing programs from one furnace to another. Even pictures can be easily and concisely managed and transferred.



Protocol

Record the performed firing cycles using the PrograBase X10 software. The recorded firing cycles can be supplemented with additional information, such as patient data, job number or material information. You can save and print the protocols or even export them as PDF files.



Download

Only one click connects PrograBase X10 with www.ivoclarvivadent.com. On the Ivoclar Vivadent website, you can download the most recent furnace software, Operating Instructions and firing tables free of charge.

System requirements

- Operating system: Windows XP[®] SP3, Windows Vista[®], Windows 7[®] or Windows 8[®] (including Microsoft .NET Framework 4*)
- Screen resolution: 1152 x 800
- Colours: min. 256
- RAM: min. 2 GB
- CPU: 1 GHz
- Min. 100 MB free memory space on the hard disk
- Internet connection

*) If Microsoft .NET Framework 4 is not yet installed on your computer, it will be installed in combination with the PrograBase software. Required memory space: approx. 300 MB.

Required accessories

USB download cable (delivered with the furnace)



This product forms part of our Fixed Prosthetics category. All the products of this category are optimally coordinated with each other.