

# Materials combination table / Massenkombinationstabelle









Chromascope shades / Chromascope Farben

	Bleach	white / weiss	yellow / gelb	light-brown / hellbraun	grey / grau	dark-brown / dunkelbraun
<b>IPS d.SIGN</b>						
<b>Opaquer</b>	 010-020    030-040	 110    120    130    140	 210    220    230    240	 310    320    330    340	 410    420    430    440	 510    520    530    540
<b>Intensiv Opaquer</b>		 white	 violet	 brown	 incisal	
<b>Margin</b>	 010-020    030-040	 110    120    130    140	 210    220    230    240	 310    320    330    340	 410    420    430    440	 510    520    530    540
<b>Intensiv Margin</b>		 yellow	 orange-pink	 add-on	 opaque	 orange
<b>Deep Dentin</b>	 010-020	 110    120    130    140	 210    220    230    240	 310    320    330    340	 410    420    430    440	 510    520    530    540
<b>Dentin</b>	 010-020    030-040	 110    120    130    140	 210    220    230    240	 310    320    330    340	 410    420    430    440	 510    520    530    540
<b>Incisal</b>	 T-S1    T-S1	 S1    S1    S1    S2	 S2    S2    S3    S3	 S3    S3    S3    S3	 S3    S3    S3    S3	 S3    S3    S3    S3
<b>Transparent</b>	 neutral	 neutral	 neutral	 neutral	 neutral	 neutral

## Materials combination table / Massenkombinationstabelle

Shades independent of any shade system / Farbsystemunabhängig

IPS d.SIGN Impulse Kit 1				
Occlusal Dentin				
	orange	brown		
Mamelon				
	light	salmon	yellow-orange	
Effect				
	Effect 1	Effect 2	Effect 3	Effect 4
Transparent				
	blue	brown-grey	orange-grey	
Incisal Edge				
	Incisal Edge			
IPS d.SIGN Impulse Kit 2				
Brilliant Dentin				
	orange	white	yellow	
Palatinal Dentin				
	red	yellow		
Insert				
	grey	orange		
Inter Incisal				
	white-blue			
Special Incisal				
	grey	yellow		
Cervical Enamel				
	khaki	orange-pink	orange	yellow

IPS d.SIGN Gingiva Kit								
Gingiva Opaquer								
	Gingiva Opaquer							
Gingiva								
	G1	G2	G3	G4	G5			
Modifier								
	GM1	GM2	GM3	GM4				
IPS d.SIGN Essence Kit								
Essence								
	brown	grey	orange	yellow				
IPS d.SIGN Stains Kit								
Stains								
	basic blue	basic red	basic yellow	black	crackliner	grey	khaki	mahogany
								
	maroon	olive	orange	vanille	white	yellow		
IPS d.SIGN Shade Kit								
Shade								
	1	2	3	4	5	6	7	